



Preparation for Cell Cycle analysis by Flow Cytometry

Reagents

- 70% ice-cold ethanol
- Trypsin
- 1X PBS
- Propidium iodide (PI)
- RNase A

Procedures

1. Harvest cells with trypsin and transfer them into a 15mL tube.
2. Sediment the cells by centrifugation at 1500 rpm for 5 min and discard the supernatant.
3. Resuspend cells in 1mL PBS and transfer the cell suspension to a 1.5 ml Eppendorf, if you have less than 1×10^6 cells.
4. Sediment the cells by centrifugation at 1500 rpm and discard the supernatant.
5. Gently resuspend the pellet in $300\mu\text{L}$ of PBS per million of cells – very important for good fixation.
6. Add slowly in drops $700\mu\text{L}$ of 100% -20°C high grade Ethanol per million of cells while slowly vortexing to prevent clumping.
7. Invert the tube a few times to insure proper mixing and a good fixation of cells.
8. Store at -20°C until day of analysis.

Cells can be stored at this point for several weeks.

9. Sediment the cells by centrifugation at 5000 rpm and remove all the ethanol.
10. Wash cells 3 times with 1X PBS.
11. Ensure all the ethanol has been removed.
12. Resuspend the pellet in 1 mL of PBS with 2.5 to $10\mu\text{g/mL}$ PI and 0.5 mg/mL RNAase A.
13. Incubate for 30 minutes at 37°C protected from light.
14. Transfer the cell suspension into a 5 mL FACS tube protected from light.